

m|devices®

AN050000 - AN059000

Suction Tubing

Sterile

INSTRUCTION FOR USE

DESCRIPTION

m|devices sterile Suction Tubing is available in smooth flexible, ribbed and special soft tubing in a variety of different lengths with one or two hexangular connectors or stepped connectors. The Suction Tubing is intended to be connected to an instrument (Yankauer, Poole, Ferguson/Frazier suction handles) to clear blood and large quantities of fluid from the surgical field during surgery, additionally can be used to remove saliva, secretions and vomit from the airway. One end of the suction tubing connects to the instrument of choice and the other is attached to a receptacle which is connected to wall suction or a stand-alone device.

Special Soft Suction Tubing is primarily used with Ferguson/Frazier suction handles when finer manipulation is required during surgery e.g. ENT and neuro procedures.

FOR USE BY A QUALIFIED CLINICIAN. THE BELOW IS ONLY A SUGGESTION AND FACILITY PROTOCOL MUST BE FOLLOWED FOR ALL CLINICAL PROCEDURES WHERE THIS PRODUCT IS USED.

CAUTION

- Single use only.
- DO NOT re-sterilise.
- DO NOT store at extreme temperatures and humidity, avoid direct sunlight. Handle with care.
- STERILE (EO), DO NOT use if the package or product has been damaged or contaminated.
- EU Notice: any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and /or patient is established.

STEPS


1. Remove product from packaging.
2. Connect the hexangular end to the suction receptacle and the other end to the selected suction handle. Some product codes have one hexangular end only to suit Ferguson/Frazier handles or to accommodate a connector to join to standard suction tubing.
3. Ensure connections are secure.
4. Test suction through the device prior to using.











m|devices

27 Llewellyn Avenue
Villawood NSW 2163
Australia
T +61(0)2 8718 2800
F +61(0)2 8718 2801

 AIRWAY
MANAGEMENT Suction Tubing

 **Multigate Co., Ltd.**
7 Lanjiang Road 312000
Shaoxing China

 **RFX+CARE International A/S**
Bakkegaardsvej 408
3050 Humlebaek Denmark

          **0123**