

m|devices®

IV012000 Secondary Infusion Set

INSTRUCTION FOR USE

DESCRIPTION

m|devices Secondary Infusion set is used to administer solution from a bag and this administration set is then added to a primary infusion line, typically to infuse small bags of antibiotics or other charted solutions. This set offers a 15micron filter vented chamber with a tubing length of 115cm.

FOR USE BY A QUALIFIED CLINICIAN. THE BELOW IS ONLY A SUGGESTION AND FACILITY PROTOCOL MUST BE FOLLOWED FOR ALL CLINICAL PROCEDURES WHERE THIS PRODUCT IS USED.

CAUTION

- Gravity use only.
- Single use only.
- DO NOT re-sterilise.
- DO NOT store at extreme temperatures and humidity, avoid direct sunlight. Handle with care.
- STERILE (EO), DO NOT use if the package or product has been damaged or contaminated.
- EU Notice: any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and /or patient is established.

STEPS

Use aseptic technique

1. Close the roller regulator.
2. Remove protective cap from spike and insert into fluid bag/bottle (If using a rigid container, open vent on spike).
3. Suspend fluid bag/bottle and squeeze lower half of the drip chamber until half full.
4. To prime line, open the roller regulator and allow fluid to flow through the tubing until all air is expelled from the line to male luer connector.
5. Close the roller regulator.
6. Remove protective cap from male luer lock connector.
7. Connect to primary infusion line by inserting male luer lock connector into access device and engaging the rotating collar.
8. Control flow rate with roller regulator.

PRECAUTION

- The security of all connections should be checked once the circuit is established and monitored during use to prevent disconnection.

m|devices

27 Llewellyn Avenue
Villawood NSW 2163
Australia
T +61(0)2 8718 2800
F +61(0)2 8718 2801

 VENOUS
ACCESS Infusion Sets

 **BQ Plus Medical Co., Ltd.**
No. 18, Che Ye Road
Che Dun Town, Songjiang
201611 Shanghai, China

 **EC REP Prolinx GmbH**
Brehmstr. 56
40239 Duesseldorf
Germany

      

     0197