

m|devices

IV020000 Transfusion Set

INSTRUCTION FOR USE

DESCRIPTION

m|devices Transfusion Sets are used to administer fluids, blood/blood products from a bag/bottle to the patient's vascular system through an intravenous catheter inserted into a vein.

This Transfusion Set consists of a tubing length of 220cm with a 200 µm filter, vented chamber and injection site.

FOR USE BY A QUALIFIED CLINICIAN. THE BELOW IS ONLY A SUGGESTION AND FACILITY PROTOCOL MUST BE FOLLOWED FOR ALL CLINICAL PROCEDURES WHERE THIS PRODUCT IS USED.

CAUTION

- Gravity use only.
- Single use only.
- DO NOT re-sterilise.
- DO NOT store at extreme temperatures and humidity, avoid direct sunlight. Handle with care.
- STERILE (EO), DO NOT use if the package or product has been damaged or contaminated.
- EU Notice: any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and /or patient is established.

STEPS

1. Close roller regulator.
2. Remove protective cap from spike and insert into fluid container (in case of rigid container, open vent on spike).
3. Suspend container.
4. Squeeze lower half of the drip chamber until half full.
5. To prime line, open the roller regulator and invert injection site until all air is expelled from line to male luer connector.
6. Close roller regulator.
7. Remove protective cap from male luer lock connector.
8. Connect and secure by inserting male luer lock connector into access device and engaging rotating collar.
9. Control flow rate with roller regulator.

USING THE INJECTION SITE


- Swab access site with antiseptic before use.
- Connect a syringe with a needleless blunt tip cannula and push directly into injection site.
- Hold the injection site whilst removing the blunt tip cannula with syringe.
- Repeat above steps each time the injection site is accessed.

PRECAUTION


- The security of all connections should be checked once the circuit is established and monitored during use to prevent disconnection.

m|devices

27 Llewellyn Avenue
Villawood NSW 2163
Australia
T +61(0)2 8718 2800
F +61(0)2 8718 2801



 **BQ Plus Medical Co., Ltd.**
No. 18, Che Ye Road
Che Dun Town, Songjiang
201611 Shanghai, China

 **Prolinx GmbH**
Brehmstr. 56
40239 Duesseldorf
Germany

