



IV020010

Gravity Transfusion Set



INSTRUCTION FOR USE

DESCRIPTION:

m|devices Gravity Transfusion Sets are used to administer fluids, blood/blood products from a bag/bottle to the patient's vascular system through an intravenous catheter inserted into a vein.

This transfusion set consists of a tubing length of 235cm with a PCA side arm, 200 µm filter, vented chamber and two swabable needleless access sites.

FOR USE BY A QUALIFIED CLINICIAN. THE BELOW IS ONLY A SUGGESTION AND FACILITY PROTOCOL MUST BE FOLLOWED FOR ALL CLINICAL PROCEDURES WHERE THIS PRODUCT IS USED.



CAUTION

DO NOT use if the package or product has been damaged or contaminated.

STEPS:

1. Close roller regulator.
2. Remove protective cap from spike and insert into fluid container (in case of rigid container, open vent on spike).
3. Suspend container.
4. Squeeze lower half of the drip chamber until half full.
5. To prime line, open the roller regulator and invert needleless access sites until all air is expelled from line including the PCA side arm to male luer connector.
6. Close roller regulator.
7. Remove protective cap from male luer lock connector.
8. Connect and secure by inserting male luer lock connector into access device and engaging rotating collar.
9. Control flow rate with roller regulator.

USING THE NEEDLELESS ACCESS SITE:

- Swab access site with antiseptic before use.
- Connect a luer slip or luer lock connector to the access site by pushing straight into the access site, whilst using a clockwise rotating motion.
- When disconnecting, hold the housing of the needleless access site to prevent disconnection from the line.
- Repeat above steps each time the needleless access site is accessed.

PRECAUTIONS:

- Do not use a needle or blunt cannula with the needleless access site.
- The security of all connections should be checked once established and monitored during use to prevent disconnection.

m|devices

27 Llewellyn Avenue
Villawood NSW 2163 Australia
T +61 (0)2 8718 2900
F +61 (0)2 8718 2801
mdevices.com.au

m|devices®