



IV072001

Closed System Transfer Set



INSTRUCTION FOR USE

DESCRIPTION

The m|devices Closed System Transfer Set is a needle free device with multi-lumen extensions designed for efficient fluid transfer.

The configuration of the set is as below:

3-Way Stopcock

3 Spikes

1 Needleless Access Site

Length 43cm

FOR USE BY A QUALIFIED CLINICIAN. THE BELOW IS ONLY A SUGGESTION AND FACILITY PROTOCOL MUST BE FOLLOWED FOR ALL CLINICAL PROCEDURES WHERE THIS PRODUCT IS USED.

CAUTION

- DO NOT use if the package or product has been damaged or contaminated.
- DO NOT store at extreme temperatures and humidity, avoid direct sunlight.
- DO NOT re-use or re-sterilise.

STEPS

1. Remove product from packaging, maintaining aseptic technique.
2. Close all click clamps and close the 3-Way tap off so that the line is closed off to the distal and side port. Remove spike caps and insert into fluid bags.
3. Open the three click clamps on the spike perforator tubing lines to allow sterile fluid to flow through the tubing to the 3-Way tap.
4. Open the click clamp to the line with the Needleless Access Site (NAS) and prime the tubing to the NAS. Once primed, close the click clamp.
5. Turn the 3-Way tap to open the line to the side port (first outlet of the 3-Way tap) to allow fluid to drip out to expel air. Once air has been removed, turn the 3-Way tap to indicate the flow towards the distal end of the set.
6. Open each of the remaining click clamps on the four distal tubing lines and allow fluid to flow through each lumen to the female luer connectors until all air has been expelled. Close click clamps.
7. If air is still visible in tubing lines, continue to prime until all air has been expelled.

PRECAUTIONS

- The security of all connections should be checked once established and monitored during use to prevent disconnection.

m|devices

27 Llewellyn Avenue
Villawood NSW 2163

Australia

T +61(0)2 8718 2800

F +61(0)2 8718 2801

m|devices[®]