



IV100001- 100003

Safety Winged Infusion Sets



INSTRUCTION FOR USE

DESCRIPTION

m|devices safety winged infusion sets shallow-angle insertion design, facilitates venepuncture of very superficial veins, e.g. hand, wrist, or scalp veins. Often selected for fragile or poorly accessible veins. The set incorporates a hypodermic needle with a needle safety shield, two bilateral flexible "wings", flexible small-bore transparent 30cm tubing with a female luer connector. The needle size options range from 19G -25G.

FOR USE BY A QUALIFIED CLINICIAN. THE BELOW IS ONLY A SUGGESTION AND FACILITY PROTOCOL MUST BE FOLLOWED FOR ALL CLINICAL PROCEDURES WHERE THIS PRODUCT IS USED.



CAUTION

DO NOT use if the package or product has been damaged or contaminated.
DO NOT store at extreme temperatures and humidity, avoid direct sunlight.

STEPS:

1. Select and prepare site and apply tourniquet as per facility protocol.
2. Using aseptic technique, open packaging and remove the sheath in a straight outward motion and inspect device.
3. Ensure the safety shield is completely pulled back to expose the needle.
4. Hold the device by the wings, pinch the wings together with thumb and fore finger to obtain a stable grip and placing your forefinger on the ribbed needle housing.
5. Anchor the vein with gentle skin traction. Insert the needle into the skin and vein at the desired angle. Blood flashback into the tubing confirms vein entry.
6. Lower the insertion angle and advance the device slightly to assure catheter entry into vein.
7. To remove the winged safety infusion set after use; Release the tourniquet and using a gauze swab apply digital pressure on the skin beyond the needle tip in the vein.
8. Slowly withdrawal needle from the vein and once completely out and away from the patient, grip the bottom of the safety shield and tubing and slide the shield upwards towards the needle tip, until the shield is fully pushed up to cover the needle and you will feel the safety shield lock in place.
9. Immediately discard locked unit into an approved sharps container/bin as per facility protocol.

PRECAUTIONS:

- Take care not to accidentally lock the safety shield in place until desired. If so, select a new winged infusion set.
- Ensure immediately after use, that the safety shield is activated/locked in place.

1800 023 321

m|devices

27 Llewellyn Avenue

Villawood NSW 2163 Australia

T +61(0)2 8718 2800

F +61(0)2 8718 2801

mdevices.com.au

m|devices®