



UR063003NS

Urine Bag Hanger

INSTRUCTION FOR USE

DESCRIPTION

m|devices Urine Bag Hanger is designed to simplify the attachment of urine bags to beds. Suitable for all major bag types.

Applying the urine bag can easily be achieved by using the two prongs on the urine hanger to attach to the corresponding perforated holes on the urine bag. The bag tubing can then be secured to the hanger using the pre-formed tubing holder that is located between the bed attachment.

CAUTION

- Single use.
- DO NOT re-sterilise.
- DO NOT store at extreme temperatures and humidity, avoid direct sunlight. Handle with care.
- NON STERILE, DO NOT use if the package or product has been damaged or contaminated.
- EU Notice: any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and /or patient is established.

STEPS

1. Attach desired drainage bag to the reinforced hanger by aligning up the two prongs on the hanger to the two pre-formed holes on the top corners of the drainage bag.
2. Secure the attachment points to ensure the bag is adequately stabilised to prevent accidental detachment.
3. Secure the urine bag tubing to the bag hanger by gently placing the tubing into pre-formed tubing holder located between the bed attachments. Check to ensure the tubing isn't kinked.
4. Attach the hanger to the bed or desired location.

PRECAUTIONS

- Ensure the drainage bag is always at or below the height of the patient's bladder.
- Check the securement of attachments prior to attaching to bed/desired location.
- Check that the urine bag tubing is kink free post placement within the pre-formed tubing holder.

m|devices

27 Llewellyn Avenue
Villawood NSW 2163
Australia
T +61 (0)2 8718 2800



 **Multigate Co., Ltd.**
4/599 Wuxie Road
Yuecheng District, Shaoxing
Zhejiang, China 312000

 **RFX+CARE International A/S**
Gammel Strandvej 22B
2990 Nivå, Denmark



Rev.04