


Standard Bore Extension Sets

Our 15cm Single Lumen Standard Bore Extension Set has a male luer lock rotating collar on one end and a NeutralSite™ on the other end.

- ⊗ NeutralSite™
- ⊗ Single Lumen
- ⊗ Male Luer lock with rotating collar
- ⊗ 3mm Tubing
- ⊗ Transparent medical grade PVC
- ⊗ Easy peel pouch



Standard Bore Extension Set				STERILE
RE-ORDER	DESCRIPTION	LENGTH	ACCESS	CARTON QTY
IV052415	Standard Bore Extension Set, Single Lumen with Male Luer Lock with Rotating Collar to NeutralSite™	15cm	NeutralSite™	200