

Standard Bore Extension Set

Single Lumen with Male Luer Lock with Rotating Collar to NeutralSite™
15cm

DESCRIPTION

The m|devices Standard Bore Extension Set is a 15cm tubing with a NeutralSite™ through to a male luer lock connection with rotating collar (RC).

FEATURES & BENEFITS

- NeutralSite™
- Standard Bore tubing
- Male luer lock end with RC
- Transparent, soft, smooth tubing
- Made from non-toxic medical grade PVC
- Peel open pouch
- Single use only


PRODUCT DATA

- Tubing diameter: 3mm
- Tubing length: 15cm
- Priming volume: 1.3mL
- Sterile EO
- Lipid resistant: No
- PSI rating: 29psi
- MRI compatible
- DEHP NO = Not detected (Result less than the Method Detection Limit EN 14372)
- Not Made with Natural Rubber Latex
- Material: PVC, PE, ABS, PC, PP, ABS, Silicone
- Shelf Life: 3 years
- ARTG: 354371

ISO STANDARDS

- EN ISO 8536-4 Infusion equipment for medical use- part 4: infusion set
- ISO 80639-7 Small-bore connectors for liquids and gases in healthcare applications — Part 7: Connectors for intravascular or hypodermic applications
- ISO 14971 Medical Devices- Application of risk management to medical device

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	NeutralSite™
2	Female Luer Lock
3	Tubing
4	Male Luer Lock
5	Cap

Standard Bore Extension Set

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV052415	Standard Bore Extension Set	15cm	200