

# S

SURGICAL  
CONSUMABLES

Irrigation  
Sets

m|devices™



We believe in effortless performance, so you can focus on caregiving.

One size never fits all. For over a decade we've been designing and manufacturing devices that simplify the roles of clinicians and carers – not just cut to a standard specification but tailored to the unique needs of Australians.

Our focus is on making performance effortless, leaving you free to focus on what matters most – caring for patients.



Leading by Design is a philosophy that reflects our commitment to research and collaboration. A crucial part of our success has been our effort to embed medical professionals in the design process. Their insight helps inform solutions that make the giving of care simple and rewarding.

The problems you face today are the solutions we're working towards for tomorrow.



# S

SURGICAL  
CONSUMABLES

## Irrigation Sets

In a surgical context you need the reliability and performance that comes from solutions carefully tailored to your needs, enabling you to perform at your best.

Our irrigation sets are designed to support an effortless maintenance of the surgical field.

## Irrigation Sets - Arthroscopy

Our Arthroscopy Irrigation Sets offer soft tubing and hand pump to simplify clearing the surgical field. They can be configured with stepped connectors with silicone adaptors or a male luer lock connector ends. Both options have a double spike configuration.

- ⊗ Length 265cm
- ⊗ Double spike Y-site configuration
- ⊗ Manual pump
- ⊗ 3 click clamps
- ⊗ Multiple connector configurations




### Irrigation Sets - Arthroscopy

RE-ORDER	DESCRIPTION	LENGTH	STERILE/NS	CARTON QTY
<b>IR010000</b>	Arthroscopy Irrigation Set	265cm	Sterile	30
<b>IR010001NS</b>	Arthroscopy Irrigation Set	265cm	Non-Sterile	50
<b>IR010002</b>	Arthroscopy Irrigation Set with MLL	265cm	Sterile	30

## Irrigation Set - Hysteroscopy

Our Hysteroscopy Irrigation Set is simple and easy to use. Each is provided with a single spike configuration, a manual hand pump made from flexible soft material with a stepped connector and a silicone cover end to connect to the hysteroscope.

- ⊗ Length 265cm
- ⊗ Single spike configuration
- ⊗ Manual pump
- ⊗ Click clamp



### Irrigation Set - Hysteroscopy

RE-ORDER	DESCRIPTION	LENGTH	STERILE
<b>IR010005</b>	Hysteroscopy Irrigation Set	265cm	30

## Irrigation Sets - Cystoscopy

Our Cystoscopy Irrigation Sets are made of high quality flexible tubing with a rigid roller regulator to accurately adjust flow rates according to clinical need during procedures.

- ⊗ Length 230cm
- ⊗ High flow rates
- ⊗ Rigid roller regulator for controlled flow
- ⊗ Soft pillow drip chamber
- ⊗ Single spike
- ⊗ Stepped connector with silicone tube
- ⊗ Available with male luer lock connector




Irrigation Sets - Cystoscopy			STERILE
RE-ORDER	DESCRIPTION	LENGTH	CARTON QTY
UR020000	Cystoscopy Irrigation Set	230cm	50
UR020002	Irrigation Set with MLL	260cm	50

## Irrigation Set - Cystoscopy/TUR

This Irrigation Set is uniquely designed for high flow irrigation post transurethral resection of the prostate or bladder to support the removal of debris, clots and tissue. The flexible soft tubing and drip chamber, double spike configuration and stepped connector with silicone tube supports the full range of clinical needs.

- ⊗ High flow rates
- ⊗ Double spike configuration
- ⊗ Soft chamber without filter
- ⊗ Two click clamps
- ⊗ Stepped connector with silicone tube



Irrigation Set - Cystoscopy/TUR			STERILE
RE-ORDER	DESCRIPTION	LENGTH	CARTON QTY
UR020001	Cystoscopy TUR Irrigation Set	240cm	50

## With you every step of the way.

As an Australian, family owned business, we're committed to supporting the best healthcare outcomes for all Australians, with Contenance Care, Venous Access and Airway Management solutions.

More than just a manufacturer of devices, we're focused on the end-to-end clinician/patient experience.

From order process, to easy disposal, we do all we can to make our role in the caregiving process seamless.

**m|devices**<sup>TM</sup>

1800 023 321



**m|devices**

27 Llewellyn Avenue  
Villawood NSW 2163  
Australia

T +61 (0)2 8718 2800

F +61 (0)2 8718 2801

[mdevices.com.au](http://mdevices.com.au)