

S

SURGICAL &
SPECIALTIES

Obstetrics

m|devices®



Designed with Midwives. Developed for Maternity Care.


At m|devices, we believe the best solutions are developed alongside the clinicians who use them every day. Created in collaboration with Australian midwives, our Cervical Dilating Balloon portfolio provides a fit-for-purpose solution for mechanical cervical ripening and labour induction.

Featuring the Cervical Dilating Balloon, Cervical Dilating Balloon Introducer and Procedure Pack, the range has been thoughtfully designed to support procedural efficiency, simplify preparation and help maternity teams focus on delivering exceptional care.



Leading by Design is a philosophy that reflects our commitment to research and collaboration. A crucial part of our success has been our effort to embed medical professionals in the design process. Their insight helps inform solutions that make the giving of care simple and rewarding.

The problems you face today are the solutions we're working towards for tomorrow.



S

SURGICAL &
SPECIALTIES

Obstetrics

In a clinical context, you need the reliability and performance that comes from solutions carefully tailored to your needs .


Our range of obstetric devices are designed to make the process of giving care simple. Their ease of use enable you to provide a reassuring patient-centred experience.

Cervical Dilating Balloon and Securement Device

Our Cervical Dilating Balloon provides a non-prescriptive alternative for midwives and obstetricians to ripen the cervix, pre-labour.

It has been designed in collaboration with clinicians and features a single balloon to improve patient comfort and the absence of eyelets creates a closed system to reduce infection risk.

- Specifically designed to ripen the cervix
- Single Balloon which can hold up to 80mL
- Avoids pharmaceutically activated hyper-stimulation
- Non-prescriptive alternative which is easy to insert
- Readily removable throughout the process
- Optional Securement Device available
- High contrast position markings
- Batch tracking sticker on packaging
- Latex free



Cervical Dilating Balloon and Securement Device			STERILE
RE-ORDER	DESCRIPTION	VOLUME	CARTON QTY
OB010020	Cervical Dilating Balloon	80mL	10
UR060001	Securement Device (8.5cm x 12.5cm)		50


Cervical Dilating Balloon Introducer

The Cervical Dilating Balloon Introducer has been designed in collaboration with clinicians to assist with the insertion of the m|devices Cervical Dilating Balloon.

Its malleable nylon material enables a simple and efficient insertion and the J-loop handle allows the placement of the Cervical Dilating Balloon with minimal movement.

The Introducer can be removed by simply reverse-bending the J-loop handle.

- Bespoke design to aid the insertion of the Cervical Dilating Balloon
- Malleable nylon material to support insertion
- J-loop handle enables the introducer to be secured in place during use
- Easily manipulated into any required angle or shape
- Designed for easy insertion and removal



Cervical Dilating Balloon Introducer			STERILE
RE-ORDER	DESCRIPTION	LENGTH	CARTON QTY
OB020000	Cervical Dilating Balloon Introducer	43.5cm	20

Cervical Dilating Procedure Pack

Efficiency

Our ‘all-in-one’ pack improves performance by eliminating collection times.

Environmental Impact

Significantly reduces individual packaging waste.

Fit for Purpose

All the necessary items are in one pack, specifically designed to support the insertion of m|devices Cervical Dilating Balloon.


- ⌋ Enables faster more efficient setup
- ⌋ Eliminates time spent collecting individual items
- ⌋ Assists with maintaining infection control principles
- ⌋ Significantly reduces individual packaging waste
- ⌋ ‘Fit for purpose’ procedure pack
- ⌋ Inclusion of the m|devices Cervical Dilating Balloon 100% Silicone, closed distal end with high contrast position markings printed on the shaft
- ⌋ Includes batch tracking sticker
- ⌋ Latex free



Cervical Dilating Balloon Procedure Pack		STERILE
RE-ORDER	DESCRIPTION	CARTON QTY
OB010002	Cervical Dilating Balloon Procedure Pack	30

Cervical Dilating Procedure Pack

Cervical Dilating Balloon Procedure Pack			STERILE		
QUANTITY	DESCRIPTION	REF	QUANTITY	DESCRIPTION	REF
1	Sterile Field 120cx120cm	A	1	m devices Cervical Dilating Balloon	G
1	Blue Rectangle Tray	B	1	XRD Gauze Swab 10x10cm 8Ply (5)	H
2	Lubricating Gel Sterile 27g	C	1	m devices Securement Device	I
4	Luer Lock Syringe 20mL	D	1	Plastic Forceps Clear	J
2	Gallipots 120mL	E	1	Tracking Label	K
1	SS Sponge Forceps 25cm	F			



When every second matters

Postpartum haemorrhage (PPH) remains a significant clinical challenge in maternity care.

- Approximately 5% of all births are affected by PPH²
- In Australia, approximately 1 in 25 vaginal births and 1 in 23 caesarean births involve a major primary PPH¹
- Many cases occur without identifiable risk factors¹

In these moments, readiness, simplicity and speed of response are essential.

Paxis™ Postpartum Balloon with Rapid Instillation Kit

Designed for critical moments in postpartum haemorrhage management.

Supporting clinicians to focus on patient care, not equipment assembly.



1. Australian Institute of Health and Welfare (AIHW), National Vore Maternity Indicators, 2003


2. NSW Health Postpartum Haemorrhage Guidelines, August 2025, GL2025_013

Paxis™ Postpartum Balloon with Rapid Instillation Kit

Paxis™ Postpartum balloon with rapid instillation kit is designed to support clinicians in the temporary control and reduction of severe postpartum uterine bleeding.

It combines a balloon catheter with an integrated rapid instillation system, helping streamline setup and enabling timely clinical response when it matters most.

- **Rapid instillation system**
Designed to support efficient fluid delivery from a sterile bag, enabling controlled balloon inflation
- **Pre-assembled configuration**
Reduces the need to gather components, supporting a smoother, more predictable setup process.
- **Soft silicone balloon (up to 500 mL)**
Designed to conform to the uterine cavity to assist in achieving haemostasis.
- **Transparent catheter**
Allows for monitoring of any ongoing bleeding



Paxis™ Postpartum Balloon with Rapid Instillation Kit			STERILE
RE-ORDER	DESCRIPTION	VOLUME	CARTON QTY
OB030000	Paxis™ Postpartum Balloon with Rapid Instillation Kit	500mL	10

Designed to support readiness

When PPH occurs, teams need to act quickly. Often in high-pressure, unpredictable situations.

Paxis™ is designed with readiness in mind:

- A pre-assembled rapid instillation system reduces setup complexity
- A single integrated kit supports efficient preparation and timely intervention

With you every step of the way.

As an Australian, family owned business, we're committed to supporting the best healthcare outcomes for all Australians, with Continence Care, Venous Access and Airway Management solutions.

More than just a manufacturer of devices, we're focused on the end-to-end clinician/patient experience.

From order process, to easy disposal, we do all we can to make our role in the caregiving process seamless.

mdevices.com.au



m|devices

27 Llewellyn Avenue
Villawood NSW 2163
Australia

T +61 (0)2 8718 2800

E info@mdevices.com.au

m|devices®

