

SaCo VLM™

Video Laryngeal Mask Airway Vision System

Product Code	Product description	Model	Size	Patient type	Patient weight (kg)	Max. ETT (mm)	Max. gastric tube (Fr)	Max .bronchoscope size (mm)	Max. OD of the instrument that can pass through the observation channel (mm)	Max OD of the instrument that can pass through airway channel	Max intracuff volume (mL)	Inner channel capacity (mL)
AN036010	SaCoVLM™ Video Larygeal Masks	MLM10C	1.0	Neonatal/infant	<5	3.5	6	2.8	3.0	4.2	10	10
AN036015	SaCoVLM™ Video Larygeal Masks	MLM15C	1.5	Neonatal/infant	5-10	4.0	8	2.8	3.0	4.2	15	13
AN036020	SaCoVLM™ Video Larygeal Masks	MLM20C	2.0	Paediatric	10-20	5.0	10	2.8	3.0	5.8	20	16
AN036025	SaCoVLM™ Video Larygeal Masks	MLM25C	2.5	Paediatric	20-30	5.5	10	2.8	3.0	5.8	25	18
AN036030	SaCoVLM™ Video Larygeal Masks	MLM30C	3.0	Adult	30-50	7.0	16	5.2	4.0	5.8	30	20
AN036040	SaCoVLM™ Video Larygeal Masks	MLM40C	4.0	Adult	50-70	7.5	16	5.5	4.0	5.8	40	23
AN036050	SaCoVLM™ Video Larygeal Masks	MLM50C	5.0	Adult	70-100	8.0	16	5.8	4.0	5.8	50	25
AN036060	SaCoVLM™ Video Larygeal Masks	MLM60C	6.0	Obesity	>100	8.0	16	5.8	4.0	5.8	60	28