

Paxis™

Post Partum Balloon Rapid Instillation Kit

DESCRIPTION

The Paxis™ Postpartum Balloon with Rapid Instillation Kit is a single-use device designed for the temporary control and reduction of severe postpartum uterine bleeding. The name Paxis™ is a blend of the latin word of 'peace' and the medical term 'haemostasis', signifying the product's goal of restoring peace by controlling postpartum haemorrhage.

The system consists of a soft silicone balloon catheter with maximum inflation volume of 500mL and a rapid instillation kit. This kit includes a pre-assembled tubing set with an IV bag spike and a three-way stopcock, designed for a quick and efficient syringe-driven inflation. The Paxis™ balloon's design allows it to conform to the shape of the uterus to apply pressure to the uterine walls and achieve haemostasis.

FEATURES & BENEFITS

Rapid Instillation Kit: Includes a pre-assembled tubing set with an IV bag spike, three-way stopcock, and a syringe to allow for rapid, syringe-driven inflation for the control of postpartum haemorrhage in an emergency.

Paxis™ Postpartum Balloon: A soft, single-use silicone balloon catheter that conforms to the uterine shape, providing effective, atraumatic haemostasis with a maximum inflation volume of 500 mL.

Transparent Catheter: Allows the clinician to easily monitor the colour and consistency of draining fluid, providing real-time feedback on the effectiveness of the treatment.

Streamlined System: The kit is designed to be ready-to-use, which simplifies the procedure and reduces preparation time in a high-stress clinical situation.

PRODUCT DATA

- Balloon Capacity: 500 mL
- Dimensions: 24 Fr, Length: 55 cm
- Sterile EO
- Shelf Life: 3 years
- ARTG: 510267

Materials:

- Balloon: Silicone
- Rapid Instillation Kit: PP, PC, PVC, ABS, PE and Rubber
- Not Made with Natural Rubber Latex

STANDARDS

- ISO 10555-4:2013 Intravascular catheters - Sterile and single-use catheters - Part 4: Balloon dilation catheters
- EN ISO 8536-10:2015- Infusion equipment for medical use- Part 10: Accessories for fluid lines for single use with pressure infusion equipment
- EN ISO 7886-1:2015- Sterile hypodermic syringes for single use- Part 1: Syringes for manual use

Paxis™ Postpartum Balloon with Rapid Instillation Kit

RE-ORDER	DESCRIPTION	BALLOON CAPACITY	CARTON/QTY
OB030000	Paxis™ Postpartum Balloon with Rapid Instillation Kit	500mL	20