



AIRWAY
MANAGEMENT

SaCo VLM™

Video Laryngeal Mask Vision
System



m|devices®
Exclusive Distributor





AIRWAY
MANAGEMENT

See the Difference Precision You Can Trust

The SaCo VLM™ Video Laryngeal Mask Vision System delivers unmatched clarity and control in airway management. Combining advanced imaging with an ergonomic design, it empowers clinicians with precision and confidence in every procedure. Engineered for efficiency, it sets a new benchmark in airway care.

SaCo VLM™ Video Laryngeal Mask Airway Vision System

Our solution includes an 8" monitor, visual stylets in three different sizes and SaCo VLM™ Video Laryngeal Mask Airway, delivering continuous, full-process visibility for reliable, advanced airway management. The system features a comprehensive range of video stylet models, purpose-designed to align seamlessly with compatible LMA sizes, ensuring optimal fit, performance and versatility across patient populations.




SaCo VLM™ Video Laryngeal Mask Airway Vision System

VS-10: Compatible with VLMA sizes #1.0 and #1.5


VS-20: Compatible with VLMA sizes #2.0 and #2.5

VS-40: Compatible with VLMA sizes #3.0, #4.0, #5.0 and #6.0

Together, these components provide clinicians with consistent visual guidance, enhanced placement accuracy, a dedicated gastric channel and greater confidence throughout the airway management process.



SaCo VLM™ Video Laryngeal Mask Airway Vision System



Core Benefits

- > **Clarity when it counts:** Real-time visualisation during insertion and intraoperative monitoring for accurate placement and reduced complications.
- > **Precision in every step:** Pre-molded design and unique distal slope enable smooth, accurate endotracheal tube (ETT) guidance.
- > **Reliable performance:** Larger contact area and optimised cuff design ensures dependable sealing and minimises reflux risk.
- > **Streamlined workflow:** Intuitive setup and dedicated channels for ventilation, gastric tube placement and visualisation, all in one integrated solution.

Clinical Applications


- > Rapid ventilation for normal and difficult airways.
- > Emergency airway management.
- > Video-assisted tracheal intubation and extubation.



Key Features

- > Visible insertion process for accurate placement and reduced air leakage.
- > Real-time image monitoring to detect reflux or cuff movement.
- > Pre-molded, ergonomic design for easier insertion and placement.
- > Unique distal slope for precise ETT insertion.
- > Dedicated channels: ventilation, gastric tube placement, visualization, inflation line and cleaning line.
- > Crystal-clear visualisation on the 8" and 10" monitor delivers enhanced airway view with quick access to photo and video capture, supported by generous storage capacity. In addition, the 10" monitor offers Wi-Fi capability and split screensharing providing enhanced visual flexibility.

SaCo VLM™ Video Laryngeal Mask Airway Vision System



Laryngeal Masks Extensive range of sizes to cater for all ages and weights

SaCo VLM™ Video Laryngeal Mask Airway Vision System			
RE-ORDER	DESCRIPTION		CARTON QTY
AN011002	UE-M8S 8" Monitor		1
AN011010	UE-M10S 10" Monitor		1
RE-ORDER	DESCRIPTION		CARTON QTY
AN038000	SaCo VLM Video Stylet - Compatible with #1.0, #1.5		1
AN038001	SaCo VLM Video Stylet - Compatible with #2.0, #2.5		1
AN038002	SaCo VLM Video Stylet - Compatible with #3.0 - #6.0		1
RE-ORDER	DESCRIPTION	SIZE	CARTON QTY
AN036010	SaCo VLM™ Video Layngeal Masks	1.0	10
AN036015	SaCo VLM™ Video Layngeal Masks	1.5	10
AN036020	SaCo VLM™ Video Layngeal Masks	2.0	10
AN036025	SaCo VLM™ Video Layngeal Masks	2.5	10
AN036030	SaCo VLM™ Video Layngeal Masks	3.0	10
AN036040	SaCo VLM™ Video Layngeal Masks	4.0	10
AN036050	SaCo VLM™ Video Layngeal Masks	5.0	10
AN036060	SaCo VLM™ Video Layngeal Masks	6.0	10

Refer to our website at www.mdevices.com.au for current stock information



SCAN HERE
TO LEARN MORE



SCAN HERE
CONTACT US FOR
A DEMO



Exciting partnership announcement

m|devices is proud to announce its partnership with UE Medical, a global airway management and anaesthesia leader.

Together we will deliver high quality solutions to healthcare professionals, offering enhanced airway management solutions that support better patient outcomes and empower clinicians with greater choice and confidence.

Through this partnership, m|devices can deliver a broader range of high-quality, clinically trusted products to the Australian market.

Discover the UE Medical Video Laryngoscope UED-A and UED-D Series.





m|devices[®]
Exclusive Distributor



mdevices.com.au

m|devices

27 Llewellyn Avenue
Villawood NSW 2163
Australia

T +61 (0)2 8718 2800

mdevices.com.au

m|devices is the trading name of
Multigate Medical Devices Pty Ltd.
ABN 69 132 290 058



A Carbon
Neutral
Company

