


S

SURGICAL &
SPECIALTIES

Paxis™
Postpartum
Balloon
with Rapid
Instillation Kit

m|devices®



When every second matters

Postpartum haemorrhage (PPH) remains a significant clinical challenge in maternity care.

- Approximately 5% of all births are affected by PPH²
- In Australia, approximately 1 in 25 vaginal births and 1 in 23 caesarean births involve a major primary PPH¹
- Many cases occur without identifiable risk factors¹

In these moments, readiness, simplicity and speed of response are essential.



S
SURGICAL &
SPECIALTIES

Paxis™ Postpartum Balloon with Rapid Instillation Kit

Designed for critical moments in postpartum haemorrhage management.

Supporting clinicians to focus on patient care, not equipment assembly.

1. Australian Institute of Health and Welfare (AIHW), National Vore Maternity Indicators, 2003


2. NSW Health Postpartum Haemorrhage Guidelines, August 2025, GL2025_013

Paxis™ Postpartum Balloon with Rapid Instillation Kit

Paxis™ Postpartum balloon with rapid instillation kit is designed to support clinicians in the temporary control and reduction of severe postpartum uterine bleeding.

It combines a balloon catheter with an integrated rapid instillation system, helping streamline setup and enabling timely clinical response when it matters most.

- Rapid instillation system**
 Designed to support efficient fluid delivery from a sterile bag, enabling controlled balloon inflation
- Pre-assembled configuration**
 Reduces the need to gather components, supporting a smoother, more predictable setup process.
- Soft silicone balloon (up to 500 mL)**
 Designed to conform to the uterine cavity to assist in achieving haemostasis.
- Transparent catheter**
 Allows for monitoring of any ongoing bleeding



Paxis™ Postpartum Balloon with Rapid Instillation Kit			STERILE
RE-ORDER	DESCRIPTION	VOLUME	CARTON QTY
OB030000	Paxis™ Postpartum Balloon with Rapid Instillation Kit	500mL	10

Designed to support readiness

When PPH occurs, teams need to act quickly. Often in high-pressure, unpredictable situations.

Paxis™ is designed with readiness in mind:

- A pre-assembled rapid instillation system reduces setup complexity
- A single integrated kit supports efficient preparation and timely intervention


See Paxis™ in action

We invite you to explore how Paxis™ can support your team:

LEARN MORE



WATCH VIDEO



REQUEST DEMO



mdevices.com.au



m|devices

27 Llewellyn Avenue
Villawood NSW 2163
Australia

T +61 (0)2 8718 2800

E info@mdevices.com.au

m|devices®

m|devices is the trading name of
Multigate Medical Devices Pty Ltd.
ABN 69 132 290 058

carbonneutral¹



A Carbon
Neutral
Company

