

S

SURGICAL &
SPECIALTIES

Theatre
Consumables

m|devices[®]



Supporting you in
delivering excellent patient
outcomes.

Achieving optimal surgical outcomes requires not only clinical expertise but also reliable, high-quality products. A well-equipped surgical suite enables clinicians to maintain focus on patient care, ensuring accuracy and efficiency in every procedure.

Our theatre consumable portfolio is designed to support surgical teams in their daily practice. From irrigation sets that provide a clear operative field, to sterile marking pens that enhance precision in pre-operative planning and intra-operative use. Each product is designed to meet the clinical demands of the surgical environment.



Contents

Irrigation Sets

06 | Arthroscopy Irrigation Sets

07 | Hysteroscopy Irrigation Set

08 | Cystoscopy Irrigation Sets

Marker Pens

12 | Skin Marker Pens

Insufflation Filter with Tubing Set

14 | Insufflation Filter with Tubing Set

S

SURGICAL &
SPECIALTIES

Reliability and effortless performance

Maximising patient outcomes in a surgical setting requires technical proficiency, excellent organisation, clear communication and products you can trust.

Our Theatre Consumables are specifically designed with this in mind.



Irrigation Sets

In a surgical context you need the reliability and performance that comes from solutions carefully tailored to your needs enabling you to perform at your best.

Our Irrigation Sets are designed to support an effortless maintenance of the surgical field.

Irrigation Sets - Arthroscopy

Our Arthroscopy Irrigation Sets offer soft tubing and hand pump to simplify clearing the surgical field. They can be configured with stepped connectors with silicone adaptors or a male luer lock connector ends. Both options have a double spike configuration.

- ⌕ Double spike Y-site configuration
- ⌕ Manual pump
- ⌕ 3 click clamps
- ⌕ Multiple connector configurations




Irrigation Sets - Arthroscopy			
RE-ORDER	DESCRIPTION	TOTAL LENGTH	CARTON QTY
IR010000	Arthroscopy Irrigation Set	265cm	30
IR010001NS*	Arthroscopy Irrigation Set	265cm	50
IR010002	Arthroscopy Irrigation Set with MLL	260cm	30

* Please note IR010001NS is a Non-Sterile Product

Irrigation Set - Hysteroscopy

Our Hysteroscopy Irrigation Set is simple and easy to use. Each is provided with a single spike configuration, a manual hand pump made from flexible soft material with a stepped connector and a silicone cover end to connect to the hysteroscope.

- Single spike configuration
- Manual pump
- Click clamp



Irrigation Set - Hysteroscopy			STERILE
RE-ORDER	DESCRIPTION	TOTAL LENGTH	CARTON QTY
IR010005	Hysteroscopy Irrigation Set	260cm	30

Irrigation Sets - Cystoscopy

Our Cystoscopy Irrigation Sets are made of high quality flexible tubing with a rigid roller regulator to accurately adjust flow rates according to clinical need during procedures.

- High flow rates
- Rigid roller regulator for controlled flow
- Soft pillow drip chamber
- Single spike
- Stepped connector with silicone tube
- Available with male luer lock connector




Irrigation Sets - Cystoscopy			STERILE
RE-ORDER	DESCRIPTION	TOTAL LENGTH	CARTON QTY
UR020000	Cystoscopy Irrigation Set	260cm	50
UR020002	Irrigation Set with MLL	285cm	50

Irrigation Set - Cystoscopy

This Cystoscopy TUR Irrigation Set is uniquely designed for high flow irrigation post transurethral resection of the prostate or bladder to support the removal of debris, clots and tissue. The flexible soft tubing and drip chamber, double spike configuration and stepped connector with silicone tube supports the full range of clinical needs.

- ⌚ High flow rates
- ⌚ Double spike configuration
- ⌚ Soft chamber without filter
- ⌚ Two click clamps
- ⌚ Stepped connector with silicone tube



Irrigation Set - Cystoscopy			STERILE
RE-ORDER	DESCRIPTION	TOTAL LENGTH	CARTON QTY
UR020001	Cystoscopy Irrigation Set	285cm	50



Where confidence meets care.

Our focus is on making performance effortless, leaving you free to focus on what matters most – caring for patients.



Skin Marking Pens

Our range of Skin Marking Pens are designed for precise marking before and during surgery, and offer continuous ink flow, bleed free marking, and fast drying.

Our focus is on making performance effortless, leaving you free to focus on surgical excellence.

Skin Marking Pens

Our fast drying, bleed free, single use marker pens are designed to deliver reliable accurate, non-toxic pre-operative marking to facilitate correct site surgery.

- Regular tip (1.0mm) and fine tip (0.5mm) options
- Scrub resistant
- Sterile and non-sterile versions
- Non-toxic
- Fast drying



Skin Marking Pens				
RE-ORDER	DESCRIPTION	LINE WIDTH	STERILE/NS	CARTON QTY
SU010111	Regular Tip Surgical Marking Pen with Ruler	1.0mm	Sterile	25
SU010111NS	Regular Tip Surgical Marking Pen with Ruler	1.0mm	Non-sterile	400
SU010121	Fine Tip Surgical Marking Pen with Ruler	0.5mm	Sterile	25



Insufflation Filters and Tubing Sets


Our range of Insufflation Filters and Tubing Sets are designed to deliver gas to patients to support intra-abdominal insufflation during laparoscopic surgeries.

Designed with exceptional filtration efficiency to protect against airborne contaminants and widely compatible with various insufflators to ease surgical set-up and use.

Irrigation Sets - Insufflation Filter with Tubing Set

m|devices Insufflation Filter with Tubing Set ensure reliable and effective carbon dioxide gas delivery to patients during minimally invasive laparoscopic surgery. These sets ensure efficient gas filtration, delivering 99.99% viral and bacterial filtration. Flat filter shape also available with luer lock connectors. All sets come with a 3m tubing set and are compatible with a variety of insufflator connectors.

- ① 99.99% filtration efficiency
- ② Withstands high gas flow
- ③ Choice of connectors
- ④ Compatible with various insufflator connectors
- ⑤ Flat shaped filter available with Luer lock connectors



Insufflation Filter with Tubing Set			STERILE
RE-ORDER	DESCRIPTION	LENGTH	CARTON QTY
SU102500	Insufflator Filter and Tubing Set with Male Luer Connector, 3m smooth	3m	25
SU102501	Insufflator Filter and Tubing Set with Push Fit Connector, 3m smooth	3m	25
SU102502	Insufflator Filter and Tubing Set with O-ring Connector, 3m smooth	3m	25
SU102503	Insufflator Filter and Tubing Set with 15/22mm Connector, 3m smooth	3m	25
Flat Filter		LENGTH	CARTON QTY
SU103000	Insufflator Flat Filter and Tubing Set with Male Luer Lock Connector, 3m smooth	3m	25
SU103100	Insufflator Flat Filter and Tubing Set with Male Luer Lock Connector, 3m ribbed	3m	25

With you every step of the way.

As an Australian, family owned business, we're committed to supporting the best healthcare outcomes for all Australians, with Continence Care, Venous Access and Airway Management solutions.

More than just a manufacturer of devices, we're focused on the end-to-end clinician/patient experience.

From order process, to easy disposal, we do all we can to make our role in the caregiving process seamless.

m|devices[®]



mdevices.com.au

m|devices

27 Llewellyn Avenue
Villawood NSW 2163
Australia

T +61 (0)2 8718 2800

mdevices.com.au

m|devices is the trading name of
Multigate Medical Devices Pty Ltd.
ABN 69 132 290 058



A Carbon
Neutral
Company

