

## Laryngoscopes - Yeescope™

The Yeescope™ one-piece, single use laryngoscope incorporates handle and blade eliminating handle and blade incompatibility and avoids the Laryngoscope becoming a source of cross contamination.

- ⊗ One-piece design avoids handle and blade incompatibility
- ⊗ Powerful illumination
- ⊗ Single use design avoids cross contamination
- ⊗ Simple function testing



| Laryngoscopes - Yeescope™ |                                     |      |             |            | NON-STERILE |
|---------------------------|-------------------------------------|------|-------------|------------|-------------|
| RE-ORDER                  | DESCRIPTION                         | SIZE | HANDLE SIZE | BLADE SIZE | CARTON QTY  |
| AN010000                  | YEESCOPE™ - Disposable Laryngoscope | MAC3 | 138mm       | 139mm      | 10          |
| AN010001                  | YEESCOPE™ - Disposable Laryngoscope | MAC4 | 140mm       | 151mm      | 10          |