

Laryngeal Mask Airways (Reinforced)

Our smooth one-piece design allows for easy insertion and the guide markings assist with positioning. Any condensation, blood secretions and regurgitated gastric contents are easily visualised through the translucent tubing. Reinforced tubing reduces kinking and supports more complex intubation.



- ⊗ Reinforced tube for difficult intubations
- ⊗ Available in adult, paediatric and neonate versions
- ⊗ Colour coded for simple size selection
- ⊗ Clear guide markers for easy positioning
- ⊗ Soft seal cuff reduces pressure and minimises potential trauma
- ⊗ A pilot balloon allows for visual and tactile assessment of the inflated / deflated cuff
- ⊗ 100% medical grade silicone
- ⊗ Translucent tubing for clear visual monitoring

Laryngeal Mask Airways (Reinforced)

STERILE

| RE-ORDER | DESCRIPTION | SIZE | COLOUR | CARTON QTY |
|----------|--|------|--------|------------|
| AN031000 | Reinforced Silicone Disposable Laryngeal Airway Mask | 1 | Purple | 10 |
| AN031001 | Reinforced Silicone Disposable Laryngeal Airway Mask | 1.5 | White | 10 |
| AN031002 | Reinforced Silicone Disposable Laryngeal Airway Mask | 2 | Red | 10 |
| AN031003 | Reinforced Silicone Disposable Laryngeal Airway Mask | 2.5 | Orange | 10 |
| AN031004 | Reinforced Silicone Disposable Laryngeal Airway Mask | 3 | Green | 10 |
| AN031005 | Reinforced Silicone Disposable Laryngeal Airway Mask | 4 | Yellow | 10 |
| AN031006 | Reinforced Silicone Disposable Laryngeal Airway Mask | 5 | Pink | 10 |