


Laryngeal Mask Airways (Reinforced)

Our smooth one-piece design allows for easy insertion and the guide markings assist with positioning. Any condensation, blood secretions and regurgitated gastric contents are easily visualised through the translucent tubing. Reinforced tubing reduces kinking and supports more complex intubation.



- ⊗ Reinforced tube for difficult intubations
- ⊗ Available in adult, paediatric and neonate versions
- ⊗ Soft seal cuff to reduce pressure and minimise trauma
- ⊗ Clear guide markers for easy positioning
- ⊗ Soft seal cuff reduces pressure and minimises potential trauma
- ⊗ A pilot balloon allows for visual and tactile assessment of the inflated / deflated cuff
- ⊗ 100% medical grade silicone
- ⊗ Translucent tubing for clear visual monitoring

Laryngeal Mask Airways (Reinforced)				STERILE
RE-ORDER	DESCRIPTION	SIZE	COLOUR	CARTON QTY
AN031000	Reinforced Silicone Disposable Laryngeal Airway Mask	1	Purple	10
AN031001	Reinforced Silicone Disposable Laryngeal Airway Mask	1.5	White	10
AN031002	Reinforced Silicone Disposable Laryngeal Airway Mask	2	Red	10
AN031003	Reinforced Silicone Disposable Laryngeal Airway Mask	2.5	Orange	10
AN031004	Reinforced Silicone Disposable Laryngeal Airway Mask	3	Green	10
AN031005	Reinforced Silicone Disposable Laryngeal Airway Mask	4	Yellow	10
AN031006	Reinforced Silicone Disposable Laryngeal Airway Mask	5	Pink	10