


## Suction Catheter Open Tip (Funnel Connector)

m|devices Suction Catheters with Funnel Connector are designed for easy passage through endotracheal and tracheostomy tubes, for the removal of airway secretions/debris that are unable to be spontaneously cleared. They feature an open distal end with an elliptical eye, soft clear tubing, incremental measurement markings and a stepped connector to connect easily to suction tubing. Available in various sizes from 6Fr through to 16Fr.

- ⊗ Open distal end with an elliptical eye
- ⊗ Stepped connector
- ⊗ Transparent
- ⊗ Incremental measurement markings
- ⊗ Available in a variety of sizes



Open Tip

Suction Catheter Open Tip (Funnel Connector)						STERILE
RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON QTY	
AN041000	Suction Catheter, Open Tip, Funnel Connector	380mm	6Fr	Light Green	50	
AN041001	Suction Catheter, Open Tip, Funnel Connector	500mm	8Fr	Light Blue	50	
AN041002	Suction Catheter, Open Tip, Funnel Connector	500mm	10Fr	Black	50	
AN041003	Suction Catheter, Open Tip, Funnel Connector	500mm	12Fr	White	50	
AN041004	Suction Catheter, Open Tip, Funnel Connector	500mm	14Fr	Green	50	
AN041005	Suction Catheter, Open Tip, Funnel Connector	500mm	16Fr	Orange	50	