


Suction Catheters Open Tip (Y Type Control Vent)

Designed to assist with a smooth insertion through endotracheal and tracheostomy tubes, for the removal of airway secretion/debris that are unable to be spontaneously cleared. The polished open catheter tip range offers soft clear tubing, incremental measurement markings commencing at 4cm from distal tip. Available in 100mm, 380mm and 560mm lengths to suit your clinical needs.

Medium style packaging also available to allow for easy storage in space limited areas.

- ⊗ Polished open tip with bilateral side eyelets
- ⊗ Stepped connector
- ⊗ Soft clear tubing
- ⊗ Colour coded to support product selection
- ⊗ Available in a wide range of diameters
- ⊗ Available in lengths 100mm, 380mm and 560mm
- ⊗ Inner sleeve with perforations to assist with maintaining aseptic technique
- ⊗ Incremental measurement markings



Open Tip

Suction Catheters Open Tip (Y Type Control Vent)

STERILE

RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON QTY
AN045002	Suction Catheter, Open Tip, Y Type Control Vent	100mm	8Fr	Light Blue	50
AN045003	Suction Catheter, Open Tip, Y Type Control Vent	100mm	10Fr	Black	50
AN045004	Suction Catheter, Open Tip, Y Type Control Vent	100mm	12Fr	White	50

Table continued on next page

m|devices[®]

m|devices
27 Llewellyn Ave
Villawood NSW 2163
Australia

T+61 (0)2 8718 2800
F+61 (0)2 8718 2801
mdevices.com.au

1800 023 321

Suction Catheters Open Tip (Y Type Control Vent)

Suction Catheters Open Tip (Y Type Control Vent)					STERILE
RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON QTY
AN042000	Suction Catheter, Open Tip, Y Type Control Vent	380mm	5Fr	Grey	50
AN042001	Suction Catheter, Open Tip, Y Type Control Vent	380mm	6Fr	Light Green	50
AN042008	Suction Catheter, Open Tip, Y Type Control Vent	380mm	7Fr	Ivory	50
AN042002	Suction Catheter, Open Tip, Y Type Control Vent	560mm	8Fr	Light Blue	50
AN042003	Suction Catheter, Open Tip, Y Type Control Vent	560mm	10Fr	Black	50
AN042004	Suction Catheter, Open Tip, Y Type Control Vent	560mm	12Fr	White	50
AN042005	Suction Catheter, Open Tip, Y Type Control Vent	560mm	14Fr	Green	50
AN042006	Suction Catheter, Open Tip, Y Type Control Vent	560mm	16Fr	Orange	50
AN042007	Suction Catheter, Open Tip, Y Type Control Vent	560mm	18Fr	Red	50
AN046000	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	380mm	5Fr	Grey	50
AN046001	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	380mm	6Fr	Light Green	50
AN046002	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	560mm	8Fr	Light Blue	50
AN046003	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	560mm	10Fr	Black	50
AN046004	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	560mm	12Fr	White	50
AN046005	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	560mm	14Fr	Green	50
AN046006	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	560mm	16Fr	Orange	50
AN046007	Suction Catheter, Open Tip, Y Type Control Vent, Medium Packaging	560mm	18Fr	Red	50