


## Suction Catheters Round Tip (Y Type Control Vent)

Designed to assist with a smooth insertion through endotracheal and tracheostomy tubes, for the removal of airway secretion/debris that are unable to be spontaneously cleared. This product comes in three sizes: 100mm, 380mm and 560mm. The round atraumatic tip range offers soft clear tubing, incremental measurement markings commencing at 4cm from distal tip.

- ⊗ Round atraumatic closed tip
- ⊗ Bilateral side eyelets
- ⊗ Available in a wide range of diameters
- ⊗ Available in lengths 100mm, 380mm and 560mm
- ⊗ Inner sleeve with perforations to assist with maintaining aseptic technique
- ⊗ Colour coded for ease of sizing selection
- ⊗ Stepped connector
- ⊗ Transparent
- ⊗ Incremental measurement markings



Round Tip

### Suction Catheters Round Tip (Y Type Control Vent)

STERILE

RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON QTY
AN044002	Suction Catheter, Round Tip, Y Type Control Vent	100mm	8Fr	Light Blue	50
AN044003	Suction Catheter, Round Tip, Y Type Control Vent	100mm	10Fr	Black	50
AN044004	Suction Catheter, Round Tip, Y Type Control Vent	100mm	12Fr	White	50

Table continued on next page

**m|devices**<sup>®</sup>

**m|devices**  
27 Llewellyn Ave  
Villawood NSW 2163  
Australia

T+61 (0)2 8718 2800  
F+61 (0)2 8718 2801  
mdevices.com.au

1800 023 321

## Suction Catheters Round Tip (Y Type Control Vent)

Suction Catheters Round Tip (Y Type Control Vent)					STERILE
RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON QTY
AN043000	Suction Catheter, Round Tip, Y Type Control Vent	380mm	5Fr	Grey	50
AN043001	Suction Catheter, Round Tip, Y Type Control Vent	380mm	6Fr	Light Green	50
AN043002	Suction Catheter, Round Tip, Y Type Control Vent	560mm	8Fr	Light Blue	50
AN043003	Suction Catheter, Round Tip, Y Type Control Vent	560mm	10Fr	Black	50
AN043004	Suction Catheter, Round Tip, Y Type Control Vent	560mm	12Fr	White	50
AN043005	Suction Catheter, Round Tip, Y Type Control Vent	560mm	14Fr	Green	50
AN043006	Suction Catheter, Round Tip, Y Type Control Vent	560mm	16Fr	Orange	50
AN043007	Suction Catheter, Round Tip, Y Type Control Vent	560mm	18Fr	Red	50