



## Suction Tubing (Non-Sterile)

Our non-sterile Suction Tubing is available in flexible, ribbed and special soft versions and comes in a variety of lengths. Each tube comes with standard hexangular connectors attached.

This product is well-suited for use in non-sterile environments such as wards.

- ⊗ Flexible, ribbed and special soft tubing options
- ⊗ Hexangular connectors
- ⊗ Kink and crush resistant tubing that maintain lumen patency
- ⊗ Variable lengths

Independent testing exceeds minimum requirements for suction tubing degree of collapse as per AS 2120.3 Clauses 5.2.2



Suction Tubing			NON-STERILE
RE-ORDER	DESCRIPTION	LENGTH	CARTON QTY
AN050000YNS	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	2m	80
AN050001YNS	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	3m	65
AN050002YNS	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	4.5m	45
AN050003YNS	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	6m	35
AN050004YNS	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	8m	30

Table continued on next page

**Suction Tubing (Non-Sterile)**

Suction Tubing			NON-STERILE
RE-ORDER	DESCRIPTION	LENGTH	CARTON QTY
<b>AN051000YNS</b>	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	2m	80
<b>AN051001YNS</b>	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	3m	65
<b>AN051002YNS</b>	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	4.5m	45
<b>AN051003YNS</b>	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	6m	35
<b>AN051004YNS</b>	Suction Tubing - Flexible, ID6mm OD9mm, 2 connectors	8m	30
<b>AN056001YNS</b>	Suction Tubing - Special Soft, ID6mm OD10mm, 1 connector	3m	65
<b>AN056006YNS</b>	Suction Tubing - Special Soft, ID6mm OD10mm, 1 connector	4m	50
<b>AN057001YNS</b>	Suction Tubing - Special Soft, ID6mm OD10mm, 2 connectors	3m	50
<b>AN059000NS</b>	Suction Tubing - Special Soft with stepped connector	0.48m	150