

Capnography Masks

Our Capnography Masks are made from non-toxic medical grade PVC and feature a 10cm extension CO2 sample line with filter and 3m oxygen tubing.

- ⊗ Adjustable strap
- ⊗ Kink resistant star lumen tubing
- ⊗ Paediatric and adult versions
- ⊗ Up to three-year shelf life
- ⊗ 'All-in-one' product to reduce inventory



| Capnography Mks | | | | STERILE |
|-------------------|---|------------|------------|---------|
| RE-ORDER | DESCRIPTION | SIZE | CARTON QTY | |
| AN065000NS | CapnographyMask - Elongated Shape • 3m Tubing | Paediatric | 50 | |
| AN065001NS | Capnography Mask - Elongated Shape, 3m Tubing | Adult | 50 | |