


Capnography Cut-Down Masks

Our Capnography Cut-Down Masks are adult size and elongated shape for maximum comfort. They feature 3m star lumen oxygen tubing and a 10cm extension CO2 sample line outlet.

Several configurations are available in this range, with options including a 0.2micron hydrophobic filter and/or 3m CO2 sample line.

The unique cut-down mask has been designed for use in ophthalmic surgery and any other procedures requiring further access to the face.

- ⊗ Cut-Down design for use in surgeries that require improved access to face
- ⊗ Contoured and comfortable to wear
- ⊗ Adjustable strap
- ⊗ Kink resistant star lumen tubing
- ⊗ Elongated mask shape to suit adult face shapes



Capnography Cut-Down Masks			NON STERILE
RE-ORDER	DESCRIPTION	SIZE	CARTON QTY
AN068000NS	Capnography Cut-Down Mask	Adult	50
AN068001NS	Capnography Cut-Down Mask, 3m CO2 Sample Line	Adult	50
AN068002NS	Capnography Cut-Down Mask with filter	Adult	50
AN068003NS	Capnography Cut-Down Mask with filter and 3m CO2 Sample Line	Adult	50

Note: All masks include 10cm Extension CO2 Sample Line Outlet, 3m Clear Oxygen Tubing