


Nasal Oxygen Cannulas

Our Nasal Oxygen Cannulas are made from non-toxic medical grade PVC and are designed to simplify oxygen administration and are soft and flexible to maximise patient comfort.

- ⊗ Kink resistant star lumen
- ⊗ Simple administration
- ⊗ Soft flexible design to maximise patient comfort

The oxygen tubing increases visibility and can assist in reducing connection errors.



Nasal Oxygen Cannulas				NON-STERILE	
RE-ORDER	DESCRIPTION	SIZE	TUBE LENGTH	CARTON QTY	
AN070000NS	Nasal Oxygen Cannula - Green Tubing	Infant	2.1m	50	
AN070001NS	Nasal Oxygen Cannula - Green Tubing	Child	2.1m	50	
AN070002NS	Nasal Oxygen Cannula - Clear Tubing	Adult	2.1m	50	
AN070003NS	Nasal Oxygen Cannula - Green Tubing	Infant	3m	50	
AN070004NS	Nasal Oxygen Cannula - Green Tubing	Child	3m	50	
AN070005NS	Nasal Oxygen Cannula - Green Tubing	Adult	3m	50	