


## Intubation Stylets

Made from malleable aluminium with a medical grade PVC sheath our intubation stylets easily bend to the required shape to assist with tube placement during difficult intubations. A straight tip design assists in tracheal tube exchange.

- ⊗ Malleable aluminium and PVC sheath design ensures stylets are easily shaped to assist in difficult intubations
- ⊗ Available in 3 sizes



### Intubation Stylets

STERILE

RE-ORDER	DESCRIPTION	SIZE	FOR ETT SIZE	CARTON QTY
AN100000	Intubating Stylet	6	3.0 - 4.5mm	20
AN100001	Intubating Stylet	10	4.5mm - 7.5mm	20
AN100002	Intubating Stylet	14	8.0 - 10.0mm	20