


## Oxygen Tubing

Our Oxygen Tubing features a channelled inner lumen to resist kinks and occlusions and is made from non-toxic medical grade PVC.

It is available in four different lengths: 2.1m, 3m, 4m and 30m.

- ⊗ Soft kink resistant tubing
- ⊗ Various lengths with optional stepped connector for the 3m length



| Oxygen Tubing |                                   |        |        | NON-STERILE |
|---------------|-----------------------------------|--------|--------|-------------|
| RE-ORDER      | DESCRIPTION                       | LENGTH | Colour | CARTON QTY  |
| AN130000NS    | Oxygen Tubing                     | 2.1m   | Clear  | 50          |
| AN130001NS    | Oxygen Tubing                     | 3m     | Clear  | 50          |
| AN130004NS    | Oxygen Tubing                     | 4m     | Clear  | 50          |
| AN130030NS    | Oxygen Tubing                     | 30m    | Clear  | 4           |
| AN135001NS    | Oxygen Tubing - Stepped Connector | 3m     | Clear  | 50          |