


## Oxygen Tubing

Our oxygen tubing features a channeled inner lumen to resist kinks and occlusions and is made from non-toxic medical grade PVC.

It is available in four different lengths: 2.1m, 3m, 4m and 30m.

- ⊗ Soft kink resistant tubing
- ⊗ Various lengths, with optional stepped connector for the 3m length



Oxygen Tubing				NON-STERILE
RE-ORDER	DESCRIPTION	LENGTH	COLOUR	CARTON QTY
AN130000NS	Oxygen Tubing	2.1m	Green	50
AN130001NS	Oxygen Tubing	3m	Clear	50
AN130004NS	Oxygen Tubing	4m	Clear	50
AN130030NS	Oxygen Tubing	30m	Green	4
AN135001NS	Oxygen Tubing - Stepped Connector	3m	Green	50

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information