



## Bubble Oxygen Tubing

Our green Bubble Oxygen Tubing has incremental sections, varying in diameter along its length. This allows it to be cut to a length which best suits connections to the oxygen source.

- ⊗ Soft kink resistant tubing
- ⊗ Green colour increases visibility and assists in reducing connection errors
- ⊗ Length 30m
- ⊗ Tubing diameter varies from 3.5mm - 7.5mm



Varying tubing diameters 3.5mm - 7.5mm

### Bubble Oxygen Tubing

NON-STERILE

RE-ORDER	DESCRIPTION	LENGTH	COLOUR	CARTON QTY
AN131030NS	Bubble Oxygen Tubing	30m	Green	4

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information