


## Bubble Oxygen Tubing

Our Bubble Oxygen Tubing has incremental sections, varying in diameter along its length. This allows it to be cut to a length which best suits connections to the oxygen source.

- ⊗ Soft kink resistant tubing
- ⊗ Length 30m
- ⊗ Tubing diameter varies from 3.5mm - 7.5mm



Varying tubing diameters 3.5mm - 7.5mm

### Bubble Oxygen T ubing

NON-STERILE

RE-ORDER	DESCRIPTION	LENGTH	Colour	CARTON QTY
<b>AN130000NS</b>	Oxygen Tubing	2.1m	Clear	50

Refer to our website at [www.mdevices.com.au](http://www.mdevices.com.au) for current stock information