

Bubble Oxygen Tubing

Our green Bubble Oxygen Tubing has incremental sections, varying in diameter along its length. This allows it to be cut to a length which best suits connections to the oxygen source.

- > Soft kink resistant tubing
- > Green colour increases visibility and assists in reducing connection errors
- > Length 30m
- > Tubing diameter varies from 3.5mm - 7.5mm



Bubble Oxygen Tubing				NON-STERILE
RE-ORDER	DESCRIPTION	LENGTH	COLOUR	CARTON QTY
AN131030NS	Bubble Oxygen Tubing	30m	Green	4