

## Anaesthetic Cushioned Masks

Our Anaesthetic Masks are made from a highly transparent, non-toxic soft PVC ensuring a comfortable facial cushion and an optimal air seal with the patients face. A colour coding system aids product selection and an inflation valve supports efficient and easy inflation and deflation. Children's masks are sweetly scented to promote calm and ease fear and anxiety.

- ⊗ Colour coded to aid product selection
- ⊗ Easy grip design
- ⊗ Facial cushion provides a strong seal
- ⊗ Children's masks are sweetly scented



Anaesthetic Cushioned Masks				NON-STERILE
RE-ORDER	DESCRIPTION	LENGTH	COLOUR	CARTON QTY
AN140001NS	Anaesthesia Mask - Neonate	1	Dark Blue	40
AN140002NS	Anaesthesia Mask - Infant - Scented	2	Blue	60
AN140003NS	Anaesthesia Mask - Child - Scented	3	Green	50
AN140004NS	Anaesthesia Mask - Small Adult	4	Yellow	40
AN140005NS	Anaesthesia Mask - Medium Adult	4.5	Orange	40
AN140005NS	Anaesthesia Mask - Adult	5	Pink	40
AN140006NS	Anaesthesia Mask - Large Adult	6	Purple	40