

Anaesthetic Circuits

Our circuits are collapsible and are available in two configurations, either corrugated with a rotatable wye connector or as Bain/Coaxial circuits with a T-connector and a 2L bag.

- ⤵ Two configurations
- ⤵ Available in three lengths
- ⤵ Adult and paediatric versions
- ⤵ Bag limb sold separately



| Anaesthetic Circuits | | | NON-STERILE |
|----------------------|-----------------------------------------------------------|--------|-------------|
| RE-ORDER | DESCRIPTION | LENGTH | CARTON QTY |
| AN150003NS | Bain Circuit with T Connector to 60cm Collapsible Circuit | 1.8m | 20 |
| AN150004NS | Bain Circuit with T Connector to 60cm Collapsible Circuit | 2.4m | 20 |
| AN150005NS | Bain Circuit with T Connector to 60cm Collapsible Circuit | 3m | 20 |
| AN150006NS | Corrugated Circuit with a Rotatable Wye Connector | 1.8m | 20 |
| AN150007NS | Corrugated Circuit with a Rotatable Wye Connector | 2.4m | 10 |
| AN150008NS | Corrugated Circuit with a Rotatable Wye Connector | 3m | 10 |
| AN150012NS | Paediatric Circuit with Wye Connector | — | 20 |
| AN151005NS | Bag Limb | 3m | 10 |