

Irrigation Sets - Arthroscopy

Our Arthroscopy Irrigation Sets offer soft tubing and hand pump to simplify clearing the surgical field. They can be configured with stepped connectors with silicone adaptors or a male luer lock connector ends. Both options have double spike configuration.

- ⊗ Length 265cm
- ⊗ Double spike Y-site configuration
- ⊗ Manual pump
- ⊗ 3 click clamps
- ⊗ Multiple connector configurations



Irrigation Sets - Arthroscopy				
RE-ORDER	DESCRIPTION	TOTAL LENGTH	STERILE/NS	CARTON QTY
IR010000	Arthroscopy Irrigation Set	265cm	Sterile	30
IR010001NS	Arthroscopy Irrigation Set	265cm	Non-Sterile	50
IR010002	Arthroscopy Irrigation Set with MLL	260cm	Sterile	50