

Infusion Sets

Infusion Sets are available with a 15 micron filter that allows for a dosage of 20 drops/mL, through 220cm kink resistant tubing. Each set comes with a needleless access site and all sets are made from medical grade PVC.

- ⊗ Intravenous solution set
- ⊗ 15 micron chamber  
20 drops per mL
- ⊗ Vented spike
- ⊗ 3mm tubing diameter
- ⊗ Male Luer lock with rotating collar
- ⊗ 220cm tube length
- ⊗ Needleless Access Site



Infusion Sets				STERILE
RE-ORDER	DESCRIPTION	LENGTH	NEEDLELESS ACCESS SITE	CARTON QTY
IV010000	Infusion Set with 15um Filter Vented Chamber and Injection Site	220cm	Latex Free Injection Site	100
IV010003	Infusion Set with 15um Filter Vented Chamber and Needleless Access Site	220cm	Needleless Access Site	50