

Infusion Sets

m|devices Infusion Set includes a vented spike attached to a drip chamber with 15µm filter roller regulator. This range provides the option of a NeutralSite™ or injection site.

The 220cm kink resistant standard tubing with male luer lock and rotating collar is made from medical grade PVC.

- ⊗ Intravenous solution set
- ⊗ 15 micron chamber 20 drops per mL
- ⊗ Vented perforator
- ⊗ 3mm tubing diameter
- ⊗ Male Luer lock with rotating collar
- ⊗ 220cm tube length
- ⊗ Peel open pouch



Infusion Sets				STERILE
RE-ORDER	DESCRIPTION	LENGTH	NEEDLELESS ACCESS SITE	CARTON QTY
IV010000	Infusion Set with 15um Filter Vented Chamber and Injection Site	220cm	Latex Free Injection Site	100
IV010003	Infusion Set with 15um Filter Vented Chamber and NeutralSite™	220cm	NeutralSite™	50