

Infusion Sets

m|devices Infusion Sets include a vented spike attached to a drip chamber with 15µm filter roller regulator to regulate flow.

Access options include: NeutralSite™, injection site or needleless access.

The kink resistant tubing is available at 220cm or 250cm. They also feature a male luer lock with rotating collar, made from medical grade PVC.

The Needleless Access Site option also offers 2 x Non Return Valves.

- ⊗ Intravenous solution set
- ⊗ 15 micron chamber 20 drops per mL
- ⊗ Vented perforator
- ⊗ 3mm tubing diameter
- ⊗ Male Luer lock with rotating collar
- ⊗ 220cm or 250cm tube length
- ⊗ Non Return Valve option available
- ⊗ Peel open pouch



Infusion Sets				STERILE
RE-ORDER	DESCRIPTION	LENGTH	NEEDLELESS ACCESS SITE	CARTON QTY
IV010000	Infusion Set with 15um Filter Vented Chamber and Injection Site	220cm	Latex Free Injection Site	100
IV010003	Infusion Set with 15um Filter Vented Chamber and NeutralSite™	220cm	NeutralSite™	50
IV010009	Infusion Set with 15um Filter Vented Chamber and 2x Needleless Access Sites & 2 x Non Return Valves	250cm	2 x Needleless Access Sites	100