

Secondary Transfusion Sets

Our Secondary Infusion Set (or Add-a-line) is made from medical grade non-toxic PVC and offers a vented spike and a short tube length of 115cm.

- ⊗ 15 micron chamber
20 drop per mL
- ⊗ Vented spike
- ⊗ 115cm tube length
- ⊗ 3mm tubing diameter
- ⊗ Male Luer lock with rotating collar



| Secondary Transfusion Sets | | | | STERILE |
|----------------------------|------------------------------------------------------------------------------|--------|------------------------|------------|
| RE-ORDER | DESCRIPTION | LENGTH | NEEDLELESS ACCESS SITE | CARTON QTY |
| IV012000 | Secondary Infusion Set with 15um Filter Vented Chamber and No Injection Site | 115cm | - | 100 |