

Secondary Infusion Set

Our Secondary Infusion Set (or Add-a-line) is made from medical grade non-toxic PVC and offers a vented spike and a short tube length of 115cm.

- ⊗ 15 micron chamber
20 drop per mL
- ⊗ Vented spike
- ⊗ 115cm tube length
- ⊗ 3mm tubing diameter
- ⊗ Male Luer lock with rotating collar



Secondary Infusion Set				STERILE
RE-ORDER	DESCRIPTION	LENGTH	NEEDLELESS ACCESS SITE	CARTON QTY
IV012000	Secondary Infusion Set with 15um Filter Vented Chamber and No Injection Site	115cm	—	100

m|devices[®]

m|devices
27 Llewellyn Ave
Villawood NSW 2163
Australia

T+61 (0)2 8718 2800
F+61 (0)2 8718 2801
mdevices.com.au

1800 023 321