


Transfusion Sets

Our Transfusion Sets lines are made from medical grade non-toxic PVC and offer a 200 micron filter that allows a dosage of 20 drops per mL.

- ⊗ 200 micron chamber (20 drops per mL)
- ⊗ Vented spike
- ⊗ Available in a range of tube lengths
- ⊗ 3mm tubing diameter
- ⊗ Male Luer lock with rotating collar
- ⊗ Available with NeutralSite™ Needleless Access Site or an option with needle injection site



Needle Injection Site



NeutralSite™



Needleless Access Site

Transfusion Sets				STERILE
RE-ORDER	DESCRIPTION	LENGTH	NEEDLELESS ACCESS SITE	CARTON QTY
IV020000	Transfusion Set with 200um Filter Vented Chamber and Injection Site	220cm	Needle Injection Site	100
IV020001	Transfusion Set with NeutralSite™	220cm	NeutralSite™	100
IV020002	Transfusion Set with Needleless Access Site	405cm	Needleless Access Site	50
IV020003	Transfusion Set with Non-Return Valve and 2 x NeutralSite™	245cm	2 x NeutralSite™	100
IV020010	Transfusion Set with 6cm PCA Side Arm, 2 x Non-Return Valves, 2 x NeutralSite™	251cm	2 x NeutralSite™	100