


## Transfusion Pump Sets

Our Transfusion Pump Sets feature a flexible dual chamber, an integral manual pump and Needleless Access Sites.

- ⊗ 200 micron chamber (20 drops per mL)
- ⊗ Manual pump
- ⊗ Range of tube lengths
- ⊗ 3mm tubing
- ⊗ Male Luer lock with rotating collar
- ⊗ Needleless Access Sites



Needleless Access Site

Transfusion Pump Sets				STERILE
RE-ORDER	DESCRIPTION	LENGTH	NEEDLELESS ACCESS SITE	CARTON QTY
IV030001	Transfusion Pump Set with Flexible Chamber and Needleless Access Site	220cm	Needleless Access Site	50
IV030003	Transfusion Pump Set with Flexible Chamber and Needleless Access Site	455cm	Needleless Access Site	50
IV030004	Transfusion Pump Set with Flexible Chamber, 2 x Needleless Access Sites and a NRV	265cm	2 x Needleless Access Site	50