

Transfusion Pump Set

Our Transfusion Pump Set operates as a gravity blood and broader solution set. It features a double spike configuration with a flexible dual chamber, an integral manual pump and Needleless Access Site.

- ⊗ 200 micron chamber (20 drops per mL)
- ⊗ Manual pump
- ⊗ 300cm tube length
- ⊗ 3mm tubing
- ⊗ Needleless Access Site



Transfusion Pump Set

STERILE

RE-ORDER	DESCRIPTION	LENGTH	NEEDLELESS ACCESS SITE	CARTON QTY
IV030008	Transfusion Pump Set with Flexible Chamber and Needleless Access Site and Luer Break	300cm	Needleless Access Site	50