


Transfusion Pump Set (Luer Break)

Our Transfusion Pump Sets operates as a gravity blood and broader solution set. It features a double spike configuration with a flexible dual chamber, an integral manual pump and Needleless Access Site.

- ⊗ 200 micron chamber (20 drops per mL)
- ⊗ Rapid infusion hand pump
- ⊗ 3mm tubing
- ⊗ Needleless Access Site
- ⊗ Male Luer lock with rotating collar
- ⊗ Luer break with female luer
- ⊗ 3 roller regulators
- ⊗ Non-vented double spike perforators



Needleless Access Site

Transfusion Pump Sets

STERILE

RE-ORDER	DESCRIPTION	LENGTH	NEEDLELESS ACCESS SITE	CARTON QTY
IV030008	Transfusion Pump Set with Flexible Chamber and Needleless Access Site and Luer Break	300cm	Needleless Access Site	50