

Burette Sets

Our In-Line Burettes have a 150mL capacity vented chamber and comes with a safety shut off valve.

- ⊗ 150mL chamber and shut off valve
- ⊗ Vented spike
- ⊗ Needleless Access Site
- ⊗ Male Luer lock connector
- ⊗ Available with administration set



Burette Sets				STERILE
RE-ORDER	DESCRIPTION	LENGTH	NEEDLELESS ACCESS SITE	CARTON QTY
IV040001	Burette - 150mL In-Line Port NAS Safety Shut off	25cm	Swabable Needleless Access Site	40
IV040002	Burette - 150mL 60 drops mL Infusion Line	230cm	2 x Swabable Needleless Access Site	25