


Burette Sets

Our in-line Burettes have a 150mL capacity vented chamber and come with a safety shut off valve. They feature; click clamp, rigid roller regulator/s for accurate flow and our needless access sites. The peel pouch facilitates easy access.

- ⊗ 150mL chamber and shut off valve
- ⊗ Vented perforator
- ⊗ Peel pouch for easy opening
- ⊗ Latex free
- ⊗ Available with administration set



Burette Sets					STERILE
RE-ORDER	DESCRIPTION	LENGTH	ACCESS	CARTON QTY	
IV040001	Burette - 150mL In-line Port NAS Safety Shut off	25cm		40	
IV040002	Burette - 150mL 60 drops mL Infusion Line	230cm	NeutralSite™	25	