


## Microbore Extensions Sets

Our Microbore Extension Sets have an inner diameter of 1mm and are available in a range of lengths.

- ⊗ Transparent, soft, smooth tubing
- ⊗ 1mm inner diameter
- ⊗ Female Luer lock to male Luer lock with rotating collar
- ⊗ Made of non-toxic medical grade PVC



Microbore Extension Sets				STERILE
RE-ORDER	DESCRIPTION	LENGTH	ACCESS	CARTON QTY
IV050000	Microbore Extension Set with Female Luer Lock to Male Luer Lock & Rotating Collar (RC)	25cm	-	100
IV050001	Microbore Extension Set with Female Luer Lock to Male Luer Lock & Rotating Collar (RC)	75cm	-	100
IV050002	Microbore Extension Set with Female Luer Lock to Male Luer Lock & Rotating Collar (RC)	150cm	-	100
IV050003	Microbore Extension Set with Female Luer Lock to Male Luer Lock & Rotating Collar (RC)	200cm	-	100

Table continued on next page

## Microbore Extensions Sets

Microbore Extension Sets				STERILE
RE-ORDER	DESCRIPTION	LENGTH	ACCESS	CARTON QTY
IV050011	Microbore Extension Set with RC Male Luer Lock to NeutralSite™ and Standard Y-Injection Site	10cm	NeutralSite™	50
IV050015	Microbore Extension Set with NeutralSite™ to Male Luer Lock with RC & Removable Slide Clamp	15cm	NeutralSite™	200
IV053001	Microbore Extension Set with 2 x Lumen and 2 x Needleless Access Site	15cm	2 x Needleless Access Site	100
IV050030	Microbore Extension Set with Female Luer Lock to Male Luer Lock with Fixed End & Slide Clamp	30cm	-	100
IV050140	Microbore Extension Set with Female Luer Lock to Male Luer Lock with Fixed End & Slide Clamp	140cm	-	100
IV050180	Microbore Extension Set with Female Luer Lock to Male Luer Lock with Fixed End & Slide Clamp	180cm	-	100