

## Microbore Extensions Sets

Our Microbore Extension Sets have an inner diameter of 1mm and are available in a range of lengths.

- ⊗ Transparent, soft, smooth tubing
- ⊗ 1mm inner diameter
- ⊗ Female Luer lock to male Luer lock with rotating collar
- ⊗ Made of non-toxic medical grade PVC



Microbore Extension Sets				STERILE
RE-ORDER	DESCRIPTION	LENGTH	NEEDLELESS ACCESS SITE	CARTON QTY
IV050000	Microbore Extension Set with Female Luer Lock to Male Luer Lock & Rotating Collar (RC)	25cm	—	100
IV050001	Microbore Extension Set with Female Luer Lock to Male Luer Lock & Rotating Collar (RC)	75cm	—	100
IV050002	Microbore Extension Set with Female Luer Lock to Male Luer Lock & Rotating Collar (RC)	150cm	—	100

Table continued on next page

Microbore Extensions Sets

Microbore Extension Sets				STERILE
RE-ORDER	DESCRIPTION	LENGTH	NEEDLELESS ACCESS SITE	CARTON QTY
IV050011	Microbore Extension Set with RC Male Luer Lock to Needleless Access Site and Standard Injection Site	10cm	Needleless Access Site	100
IV053001	Microbore Extension Set with 2 x Lumen and 2 x Needleless Access Sites	15cm	2 x Needleless Access Site	100
IV050030	Microbore Extension Set with Female Luer Lock to Male Luer Lock with Fixed End & Slide Clamp	30cm		100
IV050140	Microbore Extension Set with Female Luer Lock to Male Luer Lock with Fixed End & Slide Clamp	140cm		100
IV050180	Microbore Extension Set with Female Luer Lock to Male Luer Lock with Fixed End & Slide Clamp	180cm		100