


Microbore Extensions Sets

The m|devices Microbore Extension Set is 15cm in length with a male luer lock and rotating collar to a swabable NeutralSite™ and a removable slide clamp.

- ⊗ Microbore 1.0mm tubing
- ⊗ Priming volume: 0.2mL
- ⊗ Maximum forward pressure of 400psi
- ⊗ NeutralSite™ needleless access site
- ⊗ Not made with latex



Microbore Extension Sets				STERILE
RE-ORDER	DESCRIPTION	LENGTH	ACCESS	CARTON QTY
IV050015	Microbore Extension Set, NeutralSite™ to MLL with RC and removable Slide Clamp	15cm	NeutralSite™	200