



Minimum Volume Extension Sets

Our Minimum Volume Extension Sets have an inner diameter of 1.5mm and are available in a range of lengths.

- ⊗ Transparent, soft, smooth tubing
- ⊗ 1.5mm inner diameter
- ⊗ Female Luer lock to male Luer lock with rotating collar
- ⊗ Made of non-toxic medical grade PVC



Needle Injection Site



NeutralSite™



Needleless Access Site

Minimum Volume Extension Sets

STERILE

RE-ORDER	DESCRIPTION	LENGTH	ACCESS	CARTON QTY
IV051002	Minimum Volume Extension Set with Female Luer Lock to Male Luer Lock and Rotating Collar (RC)	600cm	-	50
IV051003	Minimum Volume Extension Set with Injection Site Female Luer Lock to Male Luer Lock & RC	165cm	Standard Injection Site	100
IV051004	Minimum Volume Extension Set with NAS Luer Lock to Male Luer Lock and RC	165cm	Needleless Access Site	100
IV051006	Minimum Volume Extension Set with NeutralSite™, Non-Return valve, Female Luer Lock to Male Luer Lock and RC	20cm	NeutralSite™	100
IV051010	Minimum Volume Extension Set with Anti-Siphon Valve	150cm	-	100